



Infrared Ear Thermometer with Probe Covers

Fast, 1 second reading



FACT BOX

Product Code:
MD1882

Product Category:
Diagnostics.

Product Approvals:
CE marked.

Power:
CR2032 battery.

Precautions:

Keep out of reach of children.

The device must stay in a stable room temperature for 15 minutes before operating. To avoid cross contamination, please clean the probe after each use.



1 Fast, 1 Sec Reading

The Medicare Infrared Ear Thermometer with Probe Covers is a fast, accurate thermometer with a user-friendly ergonomic design. It uses the infrared heat generated by the eardrum and surrounding tissue to reflect the patients body temperature. Supplied with probe covers for hygienic purposes.

FEATURES:

- Fast, 1 second measurement.
- 20 probe covers included.
- Large LCD display.
- °C/°F switchable.
- 9 memory spaces.

DIRECTIONS:

1. Install the probe cover before taking a temperature.
2. Insert the probe into the probe cover on the connection ring. Push the connection ring until you hear it 'click'.
3. Press the measure button.
4. The result will appear on the display.

RECOMMENDED USE:

For temperature monitoring.