



Digital Pulse Oximeter

Quick, non-invasive device
to measure Pulse and SpO2



FACT BOX

Product Code:
MD1913

Product Category:
Diagnostics.

Product Approvals:
CE marked.

Accuracy:

SpO2 range is 70-100% with
an accuracy of $\pm 2\%$.

Pulse rate range 25-250 with
an accuracy of $\pm 3\text{BPM}$
(beats per min).

Accessories:

x1 Pulse Oximeter
x1 Lanyard
x1 User Manual
x2 AAA Batteries
x1 Certificate Card

Power:

2 x AAA Batteries.
Auto-off function.



Low Power
Consumption



Rotate
Display



Real-time
Readings

The Medicare LifeSense® Fingertip Pulse Oximeter measures your SpO2 levels and pulse. This means it checks how well your heart is pumping oxygen through your body. It can be used to help with the monitoring of any type of condition that can affect blood oxygen levels such as COPD, Asthma, Pneumonia, Lung Cancer, Anaemia and Congenital heart defects. The Medicare LifeSense® Digital Pulse Oximeter Pulse oximetry is a quick, non-invasive and completely painless test which gives a quick and accurate results within 10 secs.

FEATURES

- OLED display.
- Pulse rate beep sound synchronised with heartbeat.
- Pulse bar graph and SpO2 waveform display.
- Four directions display.
- Track your blood oxygen levels by inputting your readings manually into the Medicare LifeSense® App.

DIRECTIONS:

Place finger in oximeter and press button. The measurement interface will appear in 3 secs & results will show in 10 secs.

RECOMMENDED USE:

Use the ring finger, middle finger or index finger to measure pulse.

Barcode: 5099390416213

MEDICARE LIFESENSE APP
Measure, Track and Share your Results

Scan to download our free Medicare LifeSense®
App and Online Health Portal or download on



@Medicarefirstaid



@Medicarefirstaid



@Medicare1staid