



Bluetooth Pulse Oximeter

Quick, non-invasive device
to measure pulse and SpO2



FACT BOX

Product Code:
MD1940


Product Category:
Diagnostics.

Product Approvals:
CE marked.

Accuracy:
SpO2 range is 70-100% with
an accuracy of $\pm 2\%$.
Pulse rate range 25-250 with
an accuracy of ± 2 BPM
(beats per min).

Accessories:
x1 Pulse Oximeter
x1 Carry Case
x1 User Manual
x2 AAA batteries

Power:
2 x AAA Batteries.
Auto-off function.



Alarm Alert
Function



Rotatable
Screen



See & Hear
Heartbeat

The Medicare LifeSense® Bluetooth Pulse Oximeter is an easy to use, non-invasive and completely painless test which gives quick and accurate results within 10 seconds. The purpose of our Bluetooth Pulse Oximeter is to measure SpO2 levels. The Bluetooth feature allows you to link the device to our Medicare LifeSense® App.

FEATURES

- Measures oxygen saturation (SpO2) and pulse rate (PR).
- See and hear heartbeat.
- Adjustable screen brightness.
- Rotatable screen with 4 direction display.
- Alarm alert function.
- Auto power off.

DIRECTIONS:

Place finger in oximeter and press the button.
The pulse results should appear within 10 seconds.

RECOMMENDED USE:

To assist with the management of conditions such as COPD, Asthma, Pneumonia, Anaemia and Congenital Heart Defects. Can also be used by athletes to monitor oxygen levels during workouts.

Barcode: 5099390417371

MEDICARE LIFESENSE APP
Measure, Track and Share your Results

Scan to download our free Medicare LifeSense®
App and Online Health Portal or download on



@Medicarefirstaid



@Medicarefirstaid



@Medicare1staid