

## Ultimate Knee Support

Recommended for people recently recovering from or looking to prevent knee injuries



### ULTIMATE



Adjustable tension straps on top and bottom offer controlled compression for weak or injured knees.



Ideal for damaged medial or lateral ligaments, cartilage damage and rehabilitation following cruciate ligament injuries.



Lined with breathable fabric which moves moisture away from your body keeping you cool and comfortable.



Padded double-layer stabilisation around the knee helps to track the patella.



Product Code:

**PX302**

Composition:

65% AIRFLEX (Neoprene), 15% Polyester, 20% Nylon.

Recognised by



Official Partner to



The Ultimate Support Range uses AIRFLEX Technology to facilitate faster recovery thanks to its therapeutic properties. AIRFLEX is an ultimate breathable material that provides the optimal level of compression and heat retention to help improve blood circulation and reduce pain.