



Hammer Toe Support

Corrects curled toes, bunions
and supports toe fractures



FACT BOX

Product Code:

MD5651S

MD5651M

MD5651L

Product Category:

Footcare

Contents:

2 x Hammer Toe Supports

Sizes:

Small (2.9cm L x 0.8cm H)

Medium (3.2cm L x 0.9cm H)

Large (3.7cm L x 1.1cm H)

Material Composition:

Medical-grade bio-gel.

x2 HAMMER
TOE SUPPORT



BIOGEL®
FORMULA



Medicare Hammer Toe Supports provide targeted relief and correction for a range of toe conditions, including curled toes, fractures, bunions, and hammertoes.

Designed with a soft, gel cushion containing mineral oil, these supports help reduce friction and pressure caused by dropped or misaligned toes.

FEATURES:

- 2 x Hammer Toe Supports
- Unique Biogel® formation eases pressure and friction

INSTRUCTIONS FOR USE:

Wear the product as shown in the picture. Remove and hand wash with warm water, then air dry in a cool place.

CAUTION:

Please keep out of reach of children. Do not use this product if you are allergic to any of its ingredients. Do not use on broken or inflamed skin. If you have diabetes, consult your doctor before use. Discontinue use and consult a doctor if you experience a rash, hives, itching, or any irritation. Store in a dry place and avoid direct sunlight.



@Medicarefirstaid



@Medicarefirstaid



@Medicare1staid