



# T2 Digital Thermometer

Perfect for fast & accurate temperature measurement

## FACT BOX

**Product Code:**  
MD061

**Product Category:**  
Diagnostics.

**Product Approvals:**  
CE marked.

**Display:**  
Liquid crystal display.

**Measurement:**  
Celsius or Fahrenheit.

**Measurement Range:**  
32.0°C - 49.9°C or  
90.0°F - 109.9°F.

**Power:**  
x1 LR41 battery.

The Medicare LifeSense® T2 Digital Thermometer is perfect for fast, accurate temperature measurement. Suitable for oral, underarm or rectal use this thermometer is waterproof and has a 10 second fast measurement time.

## FEATURES

- Auto shut-off 10 minutes after use.
- 200 hour/12 month battery life.
- Memory recall.
- Flexible tip.
- Fever alarm.

## DIRECTIONS:

Press the On/Off button next to LCD display. A tone will sound as the screen shows **100.0°C**, followed by last recorded temperature. After showing the self-test temperature, the thermometer is now in testing mode.

## RECOMMENDED USE:

Temperature monitoring.



Waterproof



Fever Notice



Flexible Tip



Memory Recall



10 Second Temperature Measurement



@Medicarefirstaid



@Medicarefirstaid



@Medicare1staid