



LIFESENSE®

A5 Blood Pressure Monitor



Reliable, accurate and easy to use

FACT BOX

Product Code:
MD1806

Product Category:
Diagnostics.

Product Approvals:
CE marked.

Cuff size: 22-42cm.

Power:
x4 AA batteries.

Colour Display:
80mm x 58.5 mm.



Irregular Heartbeat
Detection



2 x Users
60 Memories per user



WHO Risk
Category Indicator
(World Health Organisation)*



Extra Large
Full Colour Screen

Our fully automatic Medicare LifeSense® A5 Blood Pressure Monitor is designed to be effortlessly operated from the comfort of your own home and have clinically proven accuracy. The coloured screen combined with large text and buttons, make this blood pressure monitor the ideal choice for people with poor eyesight.

FEATURES

- One touch measurement.
- Large buttons & text for high visibility.
- Equipped with silent pump.
- WHO BP classification level colour indicator.
- Automatically averages last 3 measurements.
- Measure during inflation technology.
- Track your blood pressure by inputting your readings manually into the Medicare LifeSense® App.

DIRECTIONS:

Relax before taking measurement. Firmly wrap the cuff 2-3cm above bend in arm but do not wrap not too tight. Rest your arm on a table and turn on the monitor. After the results appear on the screen, you can turn off the monitor.

RECOMMENDED USE:

Recommended to monitor blood pressure.

Barcode: 5099390416688

MEDICARE LIFESENSE APP
Measure, Track and Share your Results

Scan to download our free Medicare LifeSense®
App and Online Health Portal or download on



@Medicarefirstaid



@Medicarefirstaid



@Medicare1staid