



LIFESENSE®

# A7 Blood Pressure Monitor

Lightweight, portable and easy to use



## FACT BOX

**Product Code:**  
MD1807

**Product Category:**  
Diagnostics.

**Product Approvals:**  
CE marked.

**Cuff size:** 22-42cm.

**Power:** x4 AAA batteries.

**Display:**  
Super Large Backlit Display  
(100mm x 68 mm).



Irregular Heartbeat  
Detection



2 x Users  
60 Memories per user



WHO Risk  
Category Indicator  
(World Health Organisation)\*



Extra Large  
Screen

Our lightweight Medicare LifeSense® A7 Blood Pressure Monitor is designed to connect with the Medicare LifeSense® health status monitoring app so you can effortlessly monitor your blood pressure from the comfort of your own home. The extra-large backlit screen combined with large text and buttons, make this blood pressure monitor the ideal choice for people with poor eyesight.

## FEATURES

- Bluetooth® technology.
- Large buttons & text for high visibility.
- Equipped with silent pump.
- WHO BP classification level colour indicator.
- Automatically averages last 3 measurements.
- Measure during inflation technology.

## DIRECTIONS:

Relax before taking measurement. Firmly wrap the cuff 2-3cm above bend in arm but do not wrap not too tight. Rest your arm on a table and turn on the monitor. After the results appear on the screen, you can turn off the monitor.

## RECOMMENDED USE:

Recommended to monitor blood pressure.

**Barcode:** 5099390416695

**MEDICARE LIFESENSE APP**  
Measure, Track and Share your Results

Scan to download our free Medicare LifeSense®  
App and Online Health Portal or download on



@Medicarefirstaid



@Medicarefirstaid



@Medicare1staid