



Step Counter

Step towards a healthier you



FACT BOX

Product Code:
MD1830

Product Category:
Diagnostics.

Product Approvals:
CE marked.

Power:
CR2032 battery.

Precautions:

Keep the unit out of reach of children. Do not disassemble, repair or modify the unit. Avoid exposing the unit to strong shock. Do not force open the main unit. Do not press too hard or twist the clip. The main unit is not waterproof, do not wash it or touch it with wet hands.

Using 3D technology, the Medicare Step Counter tracks steps accurately, whether it's in your pocket, on your belt or clipped to your bag. Along with a calorie consumption chart, you can clearly know how many calories you have burned.

FEATURES

- Lightweight.
- Accurate step counting.
- 3D smart sensor technology.
- Large LCD screen.
- 3D acceleration sensor.
- Reset button.

DIRECTIONS:

1. Clip the step counter to your belt or waistband for best results.
2. Open the front case and reset the data to zero.
3. Close the front case and start walking.

RECOMMENDED USE:

For managing health and calorie monitoring.



@Medicarefirstaid



@Medicarefirstaid



@Medicare1staid