



# Medicare Adjustable Hernia Aid

Provides support for reducible inguinal hernias



## FACT BOX

### Product Code:

MD540SM, MD540MD,  
MD40LG, MD540XL

### Product Category:

Home Health.

**Composition:** 55% Polyester,  
15% Rayon, 30% Rubber.

### Product Approvals:

CE marked.

**Size:** Small, Medium, Large  
and Extra Large.

**Warning:** Contains rubber  
which may cause skin  
irritations. If rash develops or  
pain persists, discontinue use  
and seek medical attention.

Medicare Adjustable Hernia Aid provides maximum support and comfort for reducible inguinal hernia for both men and women.

The flexible front panel fits the contours of the body and can be worn under clothing. Adjustable soft knitted elastic leg straps and belt hold the support firmly in place, without restricting body movement.

The removable velcro right and left support pads provide gentle gradual pressure to help reducible inguinal hernias and allows you to custom fit the support for ultimate comfort.

## FEATURES

- Support and comfort for reducible inguinal hernias.
- Suitable for both men and women.
- Adjustable fit helps allow unrestricted body movement.

## DIRECTIONS:

Measure around the middle of your hips for your size.

## RECOMMENDED USE:

Hand wash in cool water with mild detergent. Rinse thoroughly and air dry. Do not use bleach or fabric softener.