



Latex Powdered Exam Gloves

Perfect protection barrier
for your hands



FACT BOX

Product Code:
TG02

Product Category:
Infection Control.

Product Approvals:
CE marked.

Size: Medium
(100 pieces).

Precautions:

This product contains natural rubber latex. If you are allergic to latex, it may cause allergic reactions including anaphylactic reaction. Store in a cool dry place. Avoid excessive heat (30°C). Open box should be shielded from exposure to direct sunlight or fluorescent lighting. These gloves provide a biological barrier - they are not intended for use as protection from hazardous chemicals.

Mediskin Latex Powdered Exam Gloves are made from high quality natural rubber latex. The latex examination gloves are lightly powdered, essentially to make them easier to put on. The addition of powder reduces the risk of ripping or tearing them. They are also tear-resistant and offer excellent flexibility and durability.

FEATURES

- Hypoallergenic.
- Ambidextrous.
- Lightly powdered.
- Rolled cuff.
- Excellent flexibility and durability.

DIRECTIONS:

Place your hand on the template on the back of the box to find your glove size. Round up to the nearest size.

RECOMMENDED USE:

Best for medical, dental and general-purpose applications.



@Medicarefirstaid



@Medicarefirstaid



@Medicare1staid