



Latex Powder-Free Exam Gloves

Perfect protection barrier for your hands



FACT BOX

Product Code:

TG04 - Small

TG05 - Medium

TG06 - Large

TG12 - Ex. Large

Product Category:

Infection Control.

Product Approvals:

CE marked.

Precautions:

This product contains natural rubber latex, which may cause allergic reactions including anaphylactic reaction. Store in a cool dry place. Avoid excessive heat (30°C). Open box should be shielded from exposure to direct sunlight or fluorescent lighting. These gloves provide a biological barrier - they are not intended for use as protection from hazardous chemicals.

Mediskin Latex Powder-Free Exam Gloves are made from high quality natural rubber latex. The latex examination gloves are 50% thicker than powdered gloves, making them ultra tear-resistant. They are more resilient to splitting or ripping when in use. Ideal for use in a controlled environment, these powder-free gloves provide protection, great dexterity and comfort.

FEATURES

- Micro-textured for better wet and dry grip.
- Ambidextrous.
- Non-Sterile.
- Rolled cuff.
- 50% thicker than powdered gloves.

DIRECTIONS:

Place your hand on the template on the back of the box to find your glove size. Round up to the nearest size.

RECOMMENDED USE:

Best for medical, dental, general-purpose applications and for home use.



@Medicarefirstaid



@Medicarefirstaid



@Medicare1staid